## GAIN™ medium (EN)

### GAIN medium Early Stage GAIN medium Blastocyst

Cell culture media for the *in vitro* culture of human embryos and gametes

Doc.: FP09 I79 R01 C.3
Update: 16.10.2015
GAIN medium is sterilized by sterile filtration

<b>USED ABBREVIATIONS</b>
ICSI  IntraCytoplasmatic Sperm Injection
IUI  Intra-Uterine Insemination
IVF  In Vitro Fertilization
OCC  Oocyte-Corona-Cumulus

<b>INTENDED USE</b>

GAIN medium consists of 2 types of ready to use cell culture media for use with human embryos and gametes.

The media can be used for fertilization, zygotes and embryos culture through the cleavage stages and up to the expanded blastocyst stage and for embryo transfer.

- » **GAIN medium Early Stage (GAIN ES)** can be used for the handling of oocytes (in preparation of, or during IVF/ICSI) and for cell culture from fertilization up to 4-cell stage, for the processing of semen and for embryo transfer
- » **GAIN medium Blastocyst (GAIN BL)** is designed for embryo culture from 4-cell up to expanded blastocyst stage and for embryo transfer. Alternatively, the medium can be used immediately after fertilization, up till expanded blastocyst stage.

<i>For professional use only</i>
<b>COMPOSITION</b>

GAIN media are bicarbonate-buffered balanced salt solutions with 10mg/liter gentamicin and 3.5g/liter human serum albumin.

<b>MATERIAL INCLUDED WITH THE KIT</b>	
Product code	Product description
GAIN010ES	1x 10ml GAIN medium Early Stage
GAIN010BL	1x 10ml GAIN medium Blastocyst

<b>MATERIAL NOT INCLUDED WITH THE KIT</b>

- » Well dishes
- » CO<sub>2</sub> incubator (37°C – 5% CO<sub>2</sub>)
- » LAF bench (ISO Class 5)
- » Microscope
- » Syringe (e.g. 1ml Plastipack)
- » Catheter (for embryo transfer)

<b>PRODUCT SPECIFICATIONS</b>

- » Chemical composition
- » pH @ 37°C and 5% CO<sub>2</sub>: 7.20-7.45
- » Osmolality: 270-290 mOsm/kg
- » Sterility: SAL 10<sup>-3</sup>
- » Endotoxins: < 0.25 EU/ml
- » Mouse Embryo Assay (blastocysts after 96h): ≥ 80% after 30min of exposure (zygote stage)
- » Use of Ph Eur or USP grade products if applicable
- » The certificate of analysis and MSDS are available upon request

<b>PRE-USE CHECKS</b>

- » Do not use the product if it becomes cloudy, or shows any evidence of microbial contamination.
- » Do not use the product if the seal of the container is opened or defect when the product is delivered.

<b>STORAGE INSTRUCTIONS</b>

- » Suitable for transport or short term storage at elevated temperatures (up to 5 days at 37°C).
- » Store at 2-8°C, do not freeze before use.
- » Keep away from (sun)light.
- » The products can be used safely up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8°C.
- » Do not use after expiry date.

<b>WARNINGS AND PRECAUTIONS</b>

Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens. Always work under strict hygienic conditions (e.g. LAF-bench ISO Class 5) to avoid possible contamination.

GAIN media contain the antibiotic Gentamicin Sulphate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

<b>METHOD</b>

**Medium pre-equilibration**  
Incubate the medium (in culture dishes or flasks with loosened screw tops) for at least 4 hours in a CO<sub>2</sub> incubator at 37°C and 5% CO<sub>2</sub>. Ideally the medium should be pre-incubated overnight before use.

**General notes for use**  
GAIN medium should generally be used in an incubator at 37°C and 5% CO<sub>2</sub> at normal atmospheric pressure (see *Note on optimal pH and the atmospheric pressure*). The medium works well in 6 cm dishes with six 50 µl droplets covered with 6ml light paraffin oil (FertiCult Mineral Oil, FertiPro NV, Belgium). Per droplet 4 trimmed OCCs can be cultured and inseminated with approximately 0.1 million sperms.

**Standard culture method**  
» Use GAIN medium Early Stage for the handling of oocytes (in preparation of, or during IVF/ICSI), and culture of zygotes and embryos up to 4-cell stage.  
Note: After fertilization and denudation the zygotes have to be transferred to a **fresh droplet** of GAIN medium Early Stage.  
» Use GAIN medium Blastocyst for the culture of embryos from 4-cell up to expanded blastocyst stage. Make sure that the fresh medium is pre-incubated before transfer of the 4-cell stage embryos.

**Alternative culture method**  
» Use GAIN medium Early Stage for the handling of oocytes (in preparation of, or during IVF/ICSI) and fertilization.  
» Use GAIN medium Blastocyst for the culture of zygotes to expanded blastocyst stage. Make sure that the fresh medium is pre-incubated before transfer of the embryos.

**Refreshing of culture media**  
**It is advisable (but not mandatory) to transfer the embryos to fresh pre-incubated medium at least every 48 hours in culture.** Of course shorter culture times are possible (e.g. to the 8 cell stage at day 3, or morula/compaction at day 4).

**Sperm processing**  
GAIN Early Stage medium can be used to dilute and to wash semen for centrifugation. It can also be used in combination with a density gradient (Sii-Select Plus, FertiPro NV, Belgium) and/or swim-up, following standard procedures. The inseminates can be used for IUI, IVF and ICSI. Good quality sperm cells survive well over a period of 48 hours in GAIN medium.

Note: however that GAIN medium does not contain HEPES so in order to maintain a physiological pH, CO<sub>2</sub> incubation is necessary. If a HEPES-buffered medium is preferred, we recommend to use FertiCult Flushing medium (FertiPro N.V., Belgium).

**Note on optimal pH and the atmospheric pressure**  
GAIN medium must generally be used in an incubator at 37°C and 5% CO<sub>2</sub>. At normal atmospheric pressure at sea level (101.3 kPa, range: 99 - 103 kPa) this will result in a medium with a pH of approximately 7.28. At lower atmospheric pressure or higher altitudes, the percentage CO<sub>2</sub> in the incubator must be increased.

Below is a table with estimated CO<sub>2</sub> setpoint at different air pressures. The only way to obtain certainty about the pH is to measure directly under culture conditions with a suitable well-calibrated pH meter. For optimal embryo culture the pH should be: 7.28 ± 0.05.

GAIN medium contains a low concentration of Phenol Red (approx. 0.001 mM) which is sufficient to indicate the pH. At the correct pH the colour is a pale bright red. If the medium turns orange or yellowish the pH is too low, if it turns pink or purplish the pH is too high.

Altitude (meters)	Atmospheric air pressure	CO <sub>2</sub> set point*
0 (sea level)	101.3 kPa	5.0%
0 - 400 m	100 - 96 kPa	5.0 - 5.2%
400 - 800 m	96 - 92 kPa	5.3 - 5.5%
800 - 1200 m	92 - 89 kPa	5.6 - 5.8%

\* Calculated as follows: (pressure sea level / pressure altitude) x 5%.

## GAIN™ medium (FR)

### GAIN medium Early Stage GAIN medium Blastocyst

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